

DR88M Transilluminator - Parts Checklist

This DR88M transilluminator package contains the following items:

Dark Reader DR88M transilluminator (Electrical input: 12 VDC)

AC adapter (Electrical input: 120 VAC, 60 Hz)

Amber viewing screen

Pair of viewing glasses

Additional glasses and camera filters (only if ordered as separate items)

Be sure to read the Safety & Use Instructions in the accompanying Manual


viewing glasses

Use either the viewing glasses or the amber screen to view sample fluorescence.




amber screen

To view the fluorescence, place the amber screen over the sample. Or alternatively, wear the amber glasses. The amber screen can be used as a camera filter.




AC adapter



on/off switch

insert this plug into this socket



sample goes here

model number and operating voltage